

TRITCDextran

(Tetramethyl-rhodamine isothiocyanate dextran)

TdB
CONSULTANCY

Trade name: TRITC-dextran

Chemical names: Dextran(3'6'tetramethylamino dihydroxy-3oxospiro(isobenzofuran-1(3H),(9H) xanthen)-5(or 6)-yl)carbamothioate

CAS nr. NA

Structure:




Fig. 1 Structural representation of fragment of TRITC-dextran molecule.

Synthesis

Selected dextran fractions prepared from native Dextran B512F are labeled with tetramethyl-rhodamine B by a procedure similar to that described by de Belder and Granath (1). The rhodamine moiety is attached by a stable thio-carbamoyl linkage and the labeling procedure does not lead to any depolymerization of the dextran. TRITC-dextran contains from 0.001-0.008 mol. TRITC per glucose unit and at these low levels of substitution confer minimal charges to the dextran, which is an essential requirement for permeability studies.

Properties

TRITC-dextran is supplied as a red powder which dissolves freely in water or salt solutions giving a red solution. The product also dissolves in DMSO, formamide and certain other polar organic solvents but is insoluble in lower aliphatic alcohols, acetone, chloroform, dimethylformamide.

Spectral data:

Excitation is best performed at 550nm and fluorescence measured at 572 nm (see Fig.2). Studies in our laboratories have shown that the fluorescence from a TRITC-dextran solution shows only a slight increase with decreasing pH over the range pH 3-9 (see Fig.3). Similar results have been reported earlier (2.). This is of interest when making quantitative measurements. Measurements in biological media may significantly affect the fluorescence intensity which may be enhanced or depressed.




Fig. 2 Fluorescence scan of TRITC-dextran 40 in 0.025M borate pH 9.0 (9.9 mg in 50 ml buffer) Excitation 550nm; Emission 571.5nm




Fig. 3 Dependence of fluorescence of TRITC-dextran with pH in the range 4-9.

Physical chemical properties of TRITC-dextran

The dextran molecule at molecular weights greater than 5000 Daltons behaves as a flexible and extended coil in solution. Table 1 (below) shows the molecular dimensions at various molecular weights.

Dextran (Mw)	Stoke's radius
2 million	270
500 000	147
100 000	69
70	58
40	44.5
10	23.6
Albumin	35

Table 1. Molecular dimensions of dextran expressed as Stokes Radii (Å)

Dextran and TRITC-dextran will exhibit Newtonian flow characteristics i.e. the viscosity is independent of shear rate. Studies in the range pH 4 - 10 establish that the viscosity is independent of pH. The viscosities of dextran fractions at various concentrations is shown in Fig. 4.




Fig.4. The dependence of viscosity on concentration for various dextran fractions.

Storage and stability

The stability of TRITC-dextran has not been investigated in detail but it is presumed to be similar to that of FITC-dextran (see data-file). Only at elevated pH (<9) and elevated temperatures is there a risk for hydrolysis of the thiocarbonyl linkage. For further details, the reader is referred to the data-file for FITC-dextran.

Toxicity

Their toxicity patterns follow those of the parent dextrans. Clinical dextran fractions have been employed for over 50 years as plasma volume expanders.

Dextran induced anaphylactoid reactions (DIARs) have been observed in humans after injection of clinical dextran solutions (3, 4). The reactions vary from mild skin reactions to severe shock states. The incidence of severe reactions is about 1 in 2000. TRITC-dextran is also likely to display this type of behavior but few reports of problems with experimental animal have appeared.

Applications

TRITC-dextran are primarily used for studying permeability and transport in cells and tissues (5). An added benefit is that measurements of the fluorescence provide quantitative data on transport and permeability of healthy and diseased tissues. Such studies can be performed in real time by intravital fluorescence microscopy. The technique offers high sensitivity and concentrations down to 1 µg/ml can be detected in tissue fluids. Unlike FITC-dextran, the fluorescence of TRITC-dextran is not dependent on pH in the range 4 to 9. A further important property is that TRITC-dextran does not bind to artery walls (6,7).

General Procedures

The microvasculature of the hamster cheek pouch has proved to be a useful model for studying plasma leakage, e.g. following ischemia/reperfusion, exposure to inflammatory mediators (8-10). With this technique, vascular permeability can be studied in real time and be related to other microvascular events such as leukocyte adhesion and activation. The cheek pouches are examined by intravital fluorescence microscopy using suitable filters and images are captured with a digital camera (see Fig.5). A 5% TRITC-dextran 150 solution in normal saline is administered i.v. (approx. 100mg/kg bodyweight).




Fig.5 TRITC-dextran 150 injected in a hamster cheek pouch 15 min. after histamine challenge (by kind permission of E.Svensjö).

An alternative procedure using rabbit ear chambers has been described. The regenerative titanium ear-chambers (rabbits) were used to study the blood/lymph systems in the microcirculation with fluorescent dextrans. Lymph ingrowth is seen after 4-8 weeks of implantation(11).

Permeability studies using combined fluorescence stereo-microscopy were reported by Thorball(5). Extensive studies of fixation techniques in the presence of FITC-dextrans are described which are relevant also for TRITC-dextrans. Details of the microscopy set-up(filters, illumination) may be found in this article. Mullick and coworkers have described a new method for sequential quantitation of endothelial layer permeability using TRITC-dextran

4. The perfusate contained 42μg/ml TD4(12).

Permeability studies in cells

The subcellular fluorescence of FITC- and TRITC-dextrans was studied in mouse macrophage lysosomes. A comparison of the values enabled estimates of the intracellular pH. Cultures were incubated with 1 mg/ml dextran concentrations(13).

Other studies describe proton accumulation in living cells using these fluorescent probes (14). Uptake of TRITC-dextran 10 in response to osmotic cell swelling by intestine 407 cells has been described(15). TRITC-dextran 10(5mg/ml) has also been used as a fluid phase marker(16) TRITC-dextran 4 was used in studies of barrier function in prostate cancer cell cultures(17).

Mairhofer and colleagues have studied the late endosomal localization of SLP-1 in perinuclear bodies with TRITC-dextran 4(18). For studies of trans-cellular protein delivery, a TRITC-dextran 10(0.1mg/ml) was used to examine the kinetics of internalization of trans-activator fusions(19).

Permeability studies in arteries and micro-vasculature

TRITC-dextran 4 and 70 were used to study permeability in carotid arteries in folate depleted rats -concentrations of 42μg/ml were injected(21). TRITC-dextrans have also proved useful in studies of potential therapeutic treatments(22-24).

LIST OF REFERENCES

1. A.N.de Belder and K.Granath. Preparation and properties of fluorescein labelled dextrans. *Carbohydr Res.* 1973;30:375-378.
2. M.J.Geisow, Fluorescein conjugates as indicators of sub cellular pH. A critical evaluation, *Exp Cell Res.* 1984 Jan;150(1):29-35.
3. K.G.Ljungström, H.Renck, K.Strandberg et al. Ad-verse reactions to dextran in Sweden 1970-1979. *Acta Chir Scand.* 1983;149:253-262.
4. H.Hedin, W.Richter and J.Ring. Dextran-induced anaphylactoid reactions in man: role of dextran reactive antibodies. *Int Arch Allergy appl Immun.* 1976; 52:145-159.
5. N.Thorball. FITC-dextran tracers in microcircu- latory and permeability studies using combined fluorescence Stereo Microscopy, Fluorescence Stereo microscopy and electron microscopy. *Histochemistry.* 1981;71:209-233.
6. B.A.Walsh, A.E.Mullick, C.E.Banka et al. Estra- diol acts separately on the LDL-particle and artery walls to reduce LDL-accumulation. *J Lipid Res.* 2000;41:134-141.
7. G.Gardner, C.E.Banka, K.A.Roberts et al. Modified LDL- mediated increase in endothelial-layer permea- bility are attenuated with 17β-estradiol. *Arterioscler Thromb Vasc Biol.* 1999;19:854-861
8. D.Hultström and E.Svensjö. Intravital and electron microscopy study of bradykinin induced vascular permeability changes using TRITC-dextran as a tracer. *J Pathol.* 1979;129:125-133.
9. E.Svensjö. The hamster cheek pouch as a research model in inflammation. *Chapter 30 in; David she- pro(Ed) Microvascular Research-biology and Patholo- gy.* p.195-207,2006, Elsevier Academic Press.
10. E.Svensjö, E.M.Saralva, M.T.Bozza et al. Salivary gland homogenates of *Lutzomyia longipalpis* and its vasodilatory peptide maxadilan cause plasma leakage via PAC1 receptor activation. *J Vasc Res.* 2009;46:435- 446.

11. J.Jonsson, K.E.Arforss and H.C.Hint. Studies on relationships between blood and lymphatic systems within the microcirculation. *6th Europ Conf Microcirculation*. Aalborg. 1970;214-218 (Karger, Basel 1971).
12. J.Mullick, B.A.Walsh, K.M.Reiser et al. Chronic estradiol treatment attenuates stiffening glycol-oxidation and permeability of the rat carotid arteries. *Am J Physiol Heart Circ Physiol*. 2001;281:H2204-H2210.
13. M.J.Geisow. Fluorescein conjugates as indicators of subcellular pH. *Exp Cell Res*. 1984 Jan;150:29-35.
14. S.Ohkuma and B.Poole. *Proc Natl Acad Sci*. 1978 Jul;75:3327-31.
15. T.van der Wijk, S.F.Tomassen, A.B.Houtsmuller et al. Increased vesicle recycling in response to osmotic cell swelling. Cause and consequence of hypotonicity-provoked ATP release. *J Biol Chem*. 2003 Oct; 278:40020-5.
16. M.Kauppi, A.Simonsen, B.Bremnes et al. The small GTPase Rab22 interacts with EEA1 and controls endosomal membrane trafficking. *J Cell Sci*. 2002;115:899-911.
17. G.V.Shah, A.Muralidharan, M.Gokulgandhi et al. Cadherin switching and activation of β -catenin signalling underlie proinvasive actions of calcitonin-receptor axis in prostate cancer. *J Biol Chem*. 2009;284:1018-1030.
18. M.Mairhofer, M.Steiner, U.Salzer et al. Stomatolike protein-1 interacts with stomatin and is targeted to late endosomes. *J Biol Chem*. 2009;42:29218-29.
19. A.Fittipaldi, A.Ferrari, A.Zoppe et al. Cell membrane lipid rafts mediate caveolar endocytosis in HIV-1 Tat fusion proteins. *J Biol Chem*. 2003;34:141-49.
20. W.Tong, W.Dong, C.Gao et al. Charge controlled permeability of polyelectrolyte microcapsules. *J Phys Chem B*. 2005 Jul;109:13159-65.
21. J.D.Symons, A.E.Mullick, J.L.Ensunsa et al. Hyperhomocysteinemia evoked by folate depletion; effects on coronary and carotid arterial function. *Arterioscler Thromb Vasc Biol*. 2002 May;22:772-80.
22. Ngo Bich Nga Nathalie Tran, F.Knorr, W.Cheung Mak. Gradient-dependent release of the model drug TRITC-dextran from FITC-labelled BSA hydrogel microcarriers in the hair follicles of porcine ear skin. *Eur J Pharm Biopharm*. 2016 Sep 29. Epub 2016 Dec29.
23. A.Patzeit, W.C.Mak, S.Jung et al. Do nanoparticles have a future in dermal drug delivery? *J Control Release*. 2016 Sep 15; pii: S0168-3659(16)30750-7.
24. T.J.Keane, J.Dziki, E.Sobieski et al. Restoring mucosal barrier function and modifying macrophage phenotype with an extracellular matrix hydrogel: Potential therapy for ulcerative colitis. *J Crohns Colitis*. 2016 Aug 20; pii: jcw149 [Epub ahead of print].